Package 'MRMConverteR'

December 8, 2017

Type Package	
Title Convert MRM-MS .mzML files to LC-MS style .mzML files	
Version 0.0.99	
Author Tom Wilson <pre><pre><pre><pre><pre><pre><pre><pre></pre></pre></pre></pre></pre></pre></pre></pre>	
Maintainer Tom Wilson <tpw2@aber.ac.uk></tpw2@aber.ac.uk>	
Description MRMConverteR takes mzML files containing multiple MRM-MS chromatograms and converts them to LC-MS style scan data so that they can be processed by xcms.	
License GPL (>= 3)	
Depends dplyr	
Imports mzR (>= 2.1.2), purrr, tools	
<pre>URL https://github.com/wilsontom/MRMConverteR</pre>	
BugReports https://github.com/wilsontom/MRMConverteR	
LazyData TRUE	
Suggests testthat, covr, msdata	
RoxygenNote 6.0.1	
NeedsCompilation no	
R topics documented:	
convert	2
Index	3

2 convert

convert	Convert MRM-MS files to LC-MS style files	

Description

Convert MRM-MS .mzML files which contain a series of time-intensity pairs (chromatograms) to a series of LC-MS scan style m/z-intensity pairs.

Usage

```
convert(input, outpath, return = FALSE)
```

Arguments

input a character stirng of the absolute file path of a MRM-MS .mzML file

outpath a character string of the location to save converted files to. If NULL then return

must be TRUE

return logical; if TRUE then the 'peaks' list is returned from the function. Default is FALSE.

Details

Once the data is represented as a matrix of m/z and intensity pairs for each scan; the data can be further processed by open source tools such as xcms.

Author(s)

Index

convert, 2